

Product datasheet for **RC218968**

DGKA (NM_201445) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DGKA (NM_201445) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	DGKA
Synonyms:	DAGK; DAGK1; DGK-alpha
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC218968 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCCAAGGAGAGGGGCTAATAAGCCCCAGTGATTTTGCCAGCTGCAAAAATACATGGAATACTCCA
 CAAAAAGGTCAGTGATGTCTAAAGCTCTTCGAGGATGGCGAGATGGCTAAATATGTCCAAGGAGATGC
 CATTGGGTACGAGGGATTCCAGCAATTCCTGAAAATCTATCTCGAAGTGGATAATGTTCCAGACACCTA
 AGCCTGGCACTGTTTCAATCCTTTGAGACTGGTCACTGCTTAAATGAGACAAATGTGACAAAAGATGTGG
 TGTGTCTCAATGATGTTTCTGTACTTTTCCCTTCTGGAGGGTGGTCGGCCAGAAGACAAGTTAGAATT
 CACCTTCAAGCTGTACGACACGGACAGAAATGGGATCCTGGACAGCTCAGAAGTGGACAAAATTATCCTA
 CAGATGATGCGAGTGGCTGAATACCTGGATTGGGATGTGTCTGAGCTGAGGCCGATTCTTCAGGAGATGA
 TAAAAGAGATTGACTATGATGGCAGTGGCTCTGTCTCTCAAGCTGAGTGGTCCGGGCTGGGGCCACCAC
 CGTGCCACTGCTAGTGTCTGGTCTGGAGATGACTCTGAAGGACGACGGACAGCACATGTGGAGGCC
 AAGAGGTTCCCCAGACCAGTCTACTGCAATCTGTGCGAGTCAAGCATTGGTCTTGCCAAAACAGGGACTGA
 GCTGTAACCTCTGTAAGTACACTGTTACAGACCAGTGTGCCATGAAAGCCCTGCCTTGTGAAGTCAGCAC
 CTATGCCAAGTCTCGGAAGGACATTGGTGTCCAATCACATGTGTGGTGGCAGGAGGCTGTGAGTCCGGG
 CGCTGCGACCGTGTGAGAAAAGATCCGGATCTACCACAGTCTGACCGGCTGCATTGTGTATGGTGCC
 ACCTAGAGATCCACGATGACTGCCTGCAAGCGTGGGCCATGAGTGTGACTGTGGGCTGCTCCGGGATCA
 CATCTGCCTCCATCTTCCATCTATCCCAGTGTCTGGCTCTGGACCGGATCGTAAAAATAGCAAAAACA
 AGCCAGAAGACCATGGATGATTTAAATTTGAGCACCTCTGAGGCTCTGCGGATTGACCCTGTTCTAACA
 CCCACCCACTTCTCGTCTTTGTCAATCTAAGAGTGGCGGAAGCAGGGGCAGAGGGTCTCTGGAAGTT
 CCAGTATATATTAACCCCTCGACAGGTGTTCAACCTCCTAAAGGATGGTCTGAGATAGGGCTCCGATTA
 TTCAAGGATGTTCTGATAGCCGGATTTTGGTGTGGTGGAGACGGCACAGTAGGCTGGATTCTAGAGA
 CCATTGACAAAGCTAACTTGCCAGTTTTGCCTCCTGTTGTGTGTGCCCTGGTACTGGAATGATCT
 GGCTCGATGCCTAAGATGGGGAGGAGTTATGAAGGACAGAATCTGGCAAAGATCCTCAAGGATTTAGAG
 ATGAGTAAAGTGGTACATATGGATCGATGGTCTGTGGAGGTGATACCTCAACAACTGAAGAAAAAGTG
 ACCCAGTCCCCTTCAAATCATCAATAACTACTTCTCTATTGGCGTGGATGCCTCTATTGCTCATCGATT
 CCACATCATGCGAGAGAAATATCCGGAGAAGTTCAACAGCAGAATGAAGAACAAGCTATGGTACTTCGAA
 TTTGCCACATCTGAATCCATCTTCTCAACATGCAAAAAGCTGGAGGAGTCTTTGACAGTTGAGATCTGTG
 GAAACCCTGGATCTGAGCAACCTGTCCCTAGAAGGCATCGCAGTGTAAACATCCCTAGCATGCATGG
 TGGCTCCAACCTCTGGGGTGTATACCAGGAGACCCATGGGGATATCTATGGGATCAACCAGGCCCTTAGGT
 GCTACAGCTAAAGTATCACCGACCCTGATATCCTGAAAACCTGTGTACCAGACCTAAGTGACAAGAGAC
 TGGAAAGTGGTGGGCTGGAGGGTCAATTGAGATGGGCCAAATCTATACCAAGCTCAAGAATGCTGGACG
 TCGGCTGGCCAAGTGTCTGAGATCACCTTCCACACCACAAAACCCCTCCCATGCAAATGACGGAGAA
 CCCTGGATGCAGACGCCCTGTACAATCAAGATCACCCACAAGAACCAGATGCCATGCTCATGGGCCAC
 CCCCCGCTCCACCAATTTCTTTGGCTTCTTGAGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC218968 protein sequence
Red=Cloning site Green=Tags(s)

```
MAKERGLISPSDFALQKQYMEYSTKKVSDVLKLFEDGEMAKYVQGDAIGYEGFQQFLKIYLEVDNVPRLH
SLALFQSFETGHCLNETNVTKDVVCLNDVSCYFSLLEGGRPEDKLEFTFKLYDTRNGILDSSEVDKIIL
QMMRVAEYLDWDVSELRPILQEMMKEIDYDGSVSQAEWVRAGATTVPLLVLLGLEMTLKDDGQHMWRP
KRFRPRVYCNLCESSIGLGKQGLSCNLCKYTVHDQCAMKALPCEVSTYAKSRKDIGVQSHVWVRGGCESG
RCDRCQKKIRIYHSLTGLHCVWCHLEIHDDCLQAVGHECDCGLLRDHILPPSSIYPSVLASGPDRKNSKT
SQKTMDDLNLSTSEALRIDPVPNTHPLLVFVNPKSGGKQGQRVWLKWFQYILNPRQVFNLLKDGPEIGLRL
FKDVPDSRILVCGGDGTGVWILETIDKANLPVLPVAVLPLGTGNDLARCLRWGGYEGQNLAKILKDL
MSKVVHMDRWSVEVIPQQTEEKSDPVFQIINNYFSIGVDASIAHRFHIMREKYPEKFNSRMKNLWYFE
FATSEIFSTCKKLEESLVEICGKPLDLSNLSLEGIAVLNIPSMHGGSNLWGDTRRPHGDIYGINQALG
ATAKVITDPDILKTCVPLSDKRLEVVGLEGAIEMGQIYTKLKNAGRRLAKCSEITFHTTKLPMQIDGE
PWMQTPCTIKITHKNQMPMLMGPPPRSTNFFGFLS
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6044_d03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_201445


ORF Size: 2205 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_201445.1 , NP_958853.1
RefSeq Size:	2769 bp
RefSeq ORF:	2208 bp
Locus ID:	1606
UniProt ID:	P23743
Cytogenetics:	12q13.2
Protein Families:	Druggable Genome
Protein Pathways:	Glycerolipid metabolism, Glycerophospholipid metabolism, Metabolic pathways, Phosphatidylinositol signaling system
MW:	82.6 kDa
Gene Summary:	The protein encoded by this gene belongs to the eukaryotic diacylglycerol kinase family. It acts as a modulator that competes with protein kinase C for the second messenger diacylglycerol in intracellular signaling pathways. It also plays an important role in the resynthesis of phosphatidylinositols and phosphorylating diacylglycerol to phosphatidic acid. Several transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Apr 2017]


Product images:



Circular map for RC218968



Western blot validation of overexpression lysate (Cat# [LY400535]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC222395] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified DGKA protein (Cat# [TP318968]). The protein was produced from HEK293T cells transfected with DGKA cDNA clone (Cat# RC218968) using MegaTran 2.0 (Cat# [TT210002]).