

## Product datasheet for **RC218975L4V**

### ELAVL2 (NM\_004432) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ELAVL2 (NM_004432) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ELAVL2
Synonyms:	HEL-N1; HELN1; HUB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004432
ORF Size:	1077 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218975).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004432.1</a> , <a href="#">NP_004423.1</a>
RefSeq Size:	2233 bp
RefSeq ORF:	1080 bp
Locus ID:	1993
UniProt ID:	<a href="#">Q12926</a>
Cytogenetics:	9p21.3
Domains:	RRM
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 39.3 kDa

**Gene Summary:** In humans, the ELAV like RNA binding protein gene family has four members (ELAVL1-4). ELAVL RNA binding proteins recognize AU-rich elements in the 3' UTRs of gene transcripts and thereby regulate gene expression post-transcriptionally. The protein encoded by this gene binds to several 3' UTRs, including its own and also that of FOS, ID, and POU5F1. This gene encodes ELAVL2 and, like ELAVL3 and ELAVL4, is expressed specifically in neurons and primarily localizes to the cytoplasm. This protein also forms a cytosolic complex with the normally nuclear-localized ELAVL1 protein. Alternative splicing of this gene results in multiple transcript variants encoding distinct protein isoforms. [provided by RefSeq, Jul 2020]