

Product datasheet for **RC218997L4V**

DRC1 (NM_145038) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DRC1 (NM_145038) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DRC1
Synonyms:	C2orf39; CCDC164; CILD21
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145038
ORF Size:	2220 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC218997).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145038.1
RefSeq Size:	2490 bp
RefSeq ORF:	2223 bp
Locus ID:	92749
UniProt ID:	Q96MC2
Cytogenetics:	2p23.3
MW:	87 kDa


[View online »](#)

Gene Summary:

This gene encodes a central component of the nexin-dynein complex (N-DRC), which regulates the assembly of ciliary dynein. Mutations in this gene can cause ciliary dyskinesia. [provided by RefSeq, Aug 2015]