

Product datasheet for **RC219021L3V**

FRMD3 (NM_174938) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FRMD3 (NM_174938) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FRMD3
Synonyms:	4.1O; EPB41L4O; EPB41LO; P410
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_174938
ORF Size:	1791 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219021).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_174938.4
RefSeq Size:	2529 bp
RefSeq ORF:	1794 bp
Locus ID:	257019
UniProt ID:	A2A2Y4
Cytogenetics:	9q21.32
Protein Families:	Transmembrane
MW:	68.6 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a single pass membrane protein primarily found in ovaries. A similar protein in erythrocytes helps determine the shape of red blood cells, but the function of the encoded protein has not been determined. There is some evidence that this is a tumor suppressor gene, and there is also evidence linking defects in this gene to susceptibility to diabetic nephropathy in type 1 diabetes. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]