

Product datasheet for **RC219046A1V**

Human C11orf85 (MAJIN) (NM_001037225) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human C11orf85 (MAJIN) (NM_001037225) AAV Particle

Tag: Myc-DDK

Symbol: C11orf85

Synonyms: C11orf85

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC219046 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGTTTAAACCCCTTTACCTACCCGTTTCCAGAGACGAGGTTTCTTCATGCAGGACCCAATGTGTATA
AATTCAAAATCAGATATGGGAAGAGTATCAGAGGAGAAGAGATAGAAAATAAGGAAGTCATCACCCAGGA
GCTGGAGGATTCTGTCCGCGTGGTCTTGGGAACTTGGACAATCTTCAGCCCTTTGCTACAGAACCTTC
ATTGTATTTCCCTATAAAAGCAAATGGGAGAGAGTTTCCACCTGAAATTCAAACATGGGGAAATTATCT
TGATCCCTACCATTTGTTTTACTCTATATGTGGAGATGAAATGGTTCCATGAAACCTGTCACTGG
GAAACCAATAAGTGACAGTCTCTTGGGTGGTCCAGTTGAGAAGAAAGCAGTAGGAGCTGTGATGAGG
AAACGAAAACACATGGACGAGCCAGCTCCCCAGCAGGCCAGGGCTGGACAGAATAGGGAAAGAAAAAC
CCAACAAGGATTGCAGGAGACTCTGGCCTCTGATATCACTGATGTCCAGAAACAAGATTCTGAGTGGGA
CACAGCCTGCCAGGGCAATTGTCCACCCCTGCAGCACAACTCACCTCCACCTAAGGAGCGAGCAGCCA
CCGGCTTCTTTGGGTTTC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC219046 protein sequence
Red=Cloning site Green=Tags(s)

MSLKPFYTPFPETRFLHAGPNVYKFKIRYGKSIRGEEIENKEVITQELEDVSRVVLGNLDNLQPFATEHF
IVFPYKSKWERVSHLKFHGEIILIPYFVFTLYVEMKWFHENLSPGKPISDSPLGLVPVEKKAVGAVMR
KRKHMDEPSSSRPGLDRIGKEKPNKDCRRLWPLISLMSRNKILSGDTACQELSHPCSTTHLHLRSEQP
PASLGF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_001037225
ORF Size:	648 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001037225.1
RefSeq Size:	1291 bp
RefSeq ORF:	651 bp
Locus ID:	283129
UniProt ID:	Q3KP22
Cytogenetics:	11q13.1
MW:	24.8 kDa