

Product datasheet for **RC219066**

SLC38A1 (NM_030674) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | SLC38A1 (NM_030674) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | SLC38A1 |
| Synonyms: | ATA1; NAT2; SAT1; SNAT1 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

ORF Nucleotide Sequence:

>RC219066 representing NM_030674
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGATGCATTTCAAAGTGGACTCGAATTAAGTGGTTCGAAAACATGACAGTGCCCGAGGATGATAACA
 TTAGCAATGACTCCAATGATTTACCGAAGTAGAAAATGGTCAGATAAATAGCAAGTTTATTCTGATCG
 TGAAGTAGAAGAAGTCTCACAACAGCCATTTGGAAAAAAGAAGTGTGATGAGTATATCCAGGTACA
 ACCTCCTTAGGCATGTCTGTTTTAACCTAAGCAACGCCATTATGGGCAGTGGGATTTGGGACTCGCCT
 TTGCCCTGGCAAACACTGGAATCCTACTTTTTCTGGTACTTTTGACTTCAGTGACATTGTGTCTATATA
 TTCAATAAACCTCCTATTGATCTGTTCAAAGAAAACAGGCTGCATGGTGTATGAAAAGCTGGGGAAACA
 GTCTTTGGCACCACAGGAAGTTCGTAATCTTTGGAGCCACCTCTACAGAACTGGAGCAATGCTGA
 GCTACCTCTTCATCGTAAAAATGAACTACCTCTGCCATAAAGTTTCTAATGGGAAAGGAAGACATT
 TTCAGCCTGGTACGTGGATGGCCGCTTCTGGTGGTGTAGTTACCTTTGGCATAATTCTCCCTCTGTGT
 CTCTTGAAGAAGTTAGGGTATCTTGGCTATACTAGTGGATTTTCCTTGAGCTGTATGGTTTTTTTCTAA
 TTGTGGTTATTTACAAGAAATTTCAAATCCCTGCATTGTTCCAGAGCTAAATCAACAATAAGTGCTAA
 TTCAACAAATGCTGACACGTGTACGCCAAAATATGTTACCTTCAATTCAAAGACCGTGTATGCTTTACCC
 ACCATTGCATTTGCATTTGTTGGCCACCGTCAGTCTGCCAATTTACAGTGAGCTTAAAGACCGATCAC
 AGAAAAAATGCAGATGGTTTCAAACATCTCCTTTTTCGCCATGTTTGTATGTACTTCTTGACTGCCAT
 TTTTGGCTACTTGACATTCTATGACAACGTGCAGTCCGACCTCCTCACAAATATCAGAGTAAAGATGAC
 ATTCTCATCCTGACAGTGGCTGGCTGTCATTGTTGCTGTGATCCTCACAGTCCCGGTGTTATTTTTCA
 CGGTTTCGTTTCATCTTTTGAAGTGGTAAGAAAAAAGTTTAAATTTATGTCGTCATACCGTGGTTAC
 CTGCATACTCTTGGTTGTTATCAACTGTTGGTGTATCTCATACCTCCATGAAGGATATTTTTGGAGTC
 GTAGGAGTTACATCTGCTAACATGCTTATTTTCTTCTCCTTCTCTTTATTTAAAAATCACAGACC
 AGGATGGAGATAAAGGAAGTCAAAGAAATTTGGGCTGCCCTTTTCTGGGCTGGGGGTGTTGTTCTCCTT
 GGTGAGCATTCCCTTGGTGTATGACTGGGCTGCTCATCGAGTAGTGACGAAGGCCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC219066 representing NM_030674
 Red=Cloning site Green=Tags(s)

MMHFKSGLELTELQNMTPEDDNISNDSNDFTEVENQINSKFI SDRESRRSLTNSHLEKKKCD EYIPGT
 TSLGMSVFNLSNAIMGSGILGLAFALANTGILLFLVLLTSVTLLSIYSINLLLICSKETGCMVYEKLG EQ
 VFGTTGKFVIFGATSLQNTGAMLSYLFIVKNELPSAIKFLMGKEETFSAWYVDGRVLLVIVTFGIILPLC
 LLKNLGYLGYTSGFSLSCMVFFLIVVIYKFKQIPCIPELNSTISANSTNADTCTPKYVTFNSKTVYALP
 TIAFAFVCHPSVLP IYSELKDRSQKMQMVSNI SFFAMFVYFLTAIFGYLTFYDNVQSDLLHKYQSKDD
 ILILTVRLAVIVAVILTVPVLFVTVRSSLFELAKKTKFNL CRHTVVTCILLVVINLLVIFIPSMKDFIGV
 VGVTSANMLIFILPSSLYLKITDQDGDKGTQRIWAALFLGLGVLFSLVSIPLVIYDWACSSSSDEGH

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6724_c05.zip

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



* The last codon before the Stop codon of the ORF

ACCN: NM_030674

ORF Size: 1461 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030674.4](#)

RefSeq Size: 8487 bp

RefSeq ORF: 1464 bp

Locus ID: 81539

UniProt ID: [Q9H2H9](#)

Cytogenetics: 12q13.11


Domains: Aa_trans

Protein Families: Transmembrane

MW: 54 kDa

Gene Summary: Amino acid transporters play essential roles in the uptake of nutrients, production of energy, chemical metabolism, detoxification, and neurotransmitter cycling. SLC38A1 is an important transporter of glutamine, an intermediate in the detoxification of ammonia and the production of urea. Glutamine serves as a precursor for the synaptic transmitter, glutamate (Gu et al., 2001 [PubMed 11325958]).[supplied by OMIM, Mar 2008]

Product images:



Circular map for RC219066