

## Product datasheet for **RC219074L4V**

### SEC24B (NM\_001042734) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SEC24B (NM_001042734) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SEC24B
Synonyms:	SEC24
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001042734
ORF Size:	3699 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219074).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001042734.1</a> , <a href="#">NP_001036199.1</a>
RefSeq Size:	4544 bp
RefSeq ORF:	3702 bp
Locus ID:	10427
UniProt ID:	<a href="#">O95487</a>
Cytogenetics:	4q25
MW:	133.4 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the SEC24 subfamily of the SEC23/SEC24 family, which is involved in vesicle trafficking. The encoded protein is thought to be a cargo-binding component of the COPII vesicle, and is thought to be involved in the transport of secretory proteins from the endoplasmic reticulum to the Golgi apparatus. Mutations in this gene have been associated with neural tube defects, and are thought to be a result of a disruption in interactions with the protein encoded by the VANGL planar cell polarity protein 2 (VANGL2) gene. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2015]