

Product datasheet for **RC219090A1V**

Human PPIL5 (LRR1) (NM_203467) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human PPIL5 (LRR1) (NM_203467) AAV Particle

Tag: Myc-DDK

Symbol: PPIL5

Synonyms: 4-1BBLRR; LRR-1; PPIL5

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC219090 representing NM_203467
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAAGCTACACTGTGAGGTGGAGGTGATCAGCCGGCACTTGCCCGCCTTGGGGCTTAGGAACCGGGGCA
AGGGCGTCCGAGCCGTGTTGAGCCTCTGTCAGCAGACTTCCAGGAGTCAGCCGCCGTCCGAGCCTTCCT
GCTCATCTCCACCCTGAAGGACAAGCGCGGGACCCGCTATGAGCTAAGGGAGAACATTGAGCAATTCTTC
ACCAAATTTGTAGATGAGGGGAAAGCCACTGTTTCGTTAAAGGAGCCTCCTGTGGATATCTGTCTAAGTA
AGGATTCCATATGGCTCTCATATCATTCCATTCCATCTCTGCCAAGATTGGATACCGCAAAAATTTGTG
TTTGTGGAAGATTCTGTCTGAACCTTTTCATTCAAGGAAGTACTACCATGAATCTGCATTCTGTTGCCCA
CACTGTGGTCTTAGTAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC219090 representing NM_203467
Red=Cloning site Green=Tags(s)

MKLHCEVEVISRHLPALGLRNRGKGVRAVLSLCQQTSRSQPPVRAFLLISTLKDKRGTRYELRENIEQFF
TKFVDEGKATVRLKEPPVDICLSKDSIWL SYHSIPSLPRFGYRKNLCLWKILSELFHSRNYHESAFCCP
HCGLSR

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2



[View online »](#)

ACCN:	NM_203467
ORF Size:	438 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_203467.1, NP_982292.1</u>
RefSeq Size:	1049 bp
RefSeq ORF:	441 bp
Locus ID:	122769
UniProt ID:	<u>Q96L50</u>
Cytogenetics:	14q21.3
MW:	16.8 kDa