

Product datasheet for RC219133

ENSA (NM_207047) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ENSA (NM_207047) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ENSA
Synonyms: ARPP-19e
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC219133 representing NM_207047
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGGTGGTCTTGGGTGTGATGTGTGTTATTGGTTTGTAGAGGACACGCAGGAGAAAGAAGGTATTC
TGCTGAGAGAGCTGAAGAGGCAAAGCTAAAGGCCAAATACCCAAGCCTAGGACAAAAGCCTGGAGGCTC
CGACTTCTCATGAAGAGACTCCAGAAAGGGCAAAGTACTTTGACTCAGGAGACTACAACATGGCCAAA
GCCAAGATGAAGAATAAGCAGCTGCCAAGTGACAGGACCAGACAAGAACCTGGTGACTGGTGATCACATCC
CCACCCACAGGATCTGCCCCAGAGAAAGTCTCGCTCGTACCAGCAAGCTTGGCGGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC219133 representing NM_207047
Red=Cloning site Green=Tags(s)
MAGGLGCDVCYWFVEDTQEKEGILPERAEEAKLKAKYPSLGQKPGGSDFLMKRLQKQKYFDSGDYNAK
AKMKNKQLPSAGPDKNLVTGDHIPTQDLPQRKSSLVTSKLAG
TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_207047

ORF Size: 339 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_207047.1</u> , <u>NP_996930.1</u>
RefSeq Size:	2747 bp
RefSeq ORF:	342 bp
Locus ID:	2029
UniProt ID:	<u>O43768</u>
Cytogenetics:	1q21.3
Protein Families:	Druggable Genome
MW:	12.2 kDa
Gene Summary:	<p>The protein encoded by this gene belongs to a highly conserved cAMP-regulated phosphoprotein (ARPP) family. This protein was identified as an endogenous ligand for the sulfonylurea receptor, ABCC8/SUR1. ABCC8 is the regulatory subunit of the ATP-sensitive potassium (KATP) channel, which is located on the plasma membrane of pancreatic beta cells and plays a key role in the control of insulin release from pancreatic beta cells. This protein is thought to be an endogenous regulator of KATP channels. In vitro studies have demonstrated that this protein modulates insulin secretion through the interaction with KATP channel, and this gene has been proposed as a candidate gene for type 2 diabetes. At least eight alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]</p>