

Product datasheet for RC219141L2V

VNN3 (NM_018399) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	VNN3 (NM_018399) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VNN3
Synonyms:	HSA238982; MGC124285; MGC171203; OTTMUSP00000022908; PAGEL-beta; PAGEL-eta; PAGEL-zeta; vanin 3; vascular non-inflammatory molecule 3
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_018399
ORF Size:	822 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219141).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018399.3 , NP_060869.2
RefSeq Size:	1733 bp
RefSeq ORF:	824 bp
Locus ID:	55350
Cytogenetics:	6q23.2
Protein Families:	Transmembrane
MW:	28.4 kDa


[View online »](#)

Gene Summary:

This gene is the central gene in a cluster of three vanin genes on chromosome 6q23-q24. Extensive alternative splicing has been described; the two most common variants are represented as RefSeqs. [provided by RefSeq, Apr 2014]