

Product datasheet for RC219196L4V

PLEKHB2 (NM_001100623) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PLEKHB2 (NM_001100623) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PLEKHB2
Synonyms:	EVT2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001100623
ORF Size:	666 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219196).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001100623.1 , NP_001094093.1
RefSeq Size:	4326 bp
RefSeq ORF:	669 bp
Locus ID:	55041
UniProt ID:	Q96CS7
Cytogenetics:	2q21.1
MW:	24.6 kDa
Gene Summary:	Involved in retrograde transport of recycling endosomes.[UniProtKB/Swiss-Prot Function]



[View online »](#)