

Product datasheet for RC219208L4

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

LAMP1 (NM_005561) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: LAMP1 (NM_005561) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: LAMP1

Synonyms: CD107a; LAMPA; LGP120

Mammalian Cell Puromycin

Selection:

Vector:

pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC219208).

Sgfl-Mlul

Sequence:

Restriction Sites:

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_005561

ORF Size: 1251 bp





LAMP1 (NM_005561) Human Tagged Lenti ORF Clone - RC219208L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 005561.3</u>, <u>NP 005552.3</u>

 RefSeq Size:
 2560 bp

 RefSeq ORF:
 1254 bp

 Locus ID:
 3916

 UniProt ID:
 P11279

Cytogenetics: 13q34

Domains:

Protein Families: ES Cell Differentiation/IPS, Transmembrane

Lamp

Protein Pathways: Lysosome MW: 44.88 kDa


Gene Summary: The protein encoded by this gene is a member of a family of membrane glycoproteins. This

glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell


metastasis. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC219208L4



Double digestion of RC219208L4 using Sgfl and Mlul $\,$