

## Product datasheet for **RC219227L3V**

### PRPS2 (NM\_001039091) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PRPS2 (NM_001039091) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PRPS2
Synonyms:	PRSII
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001039091
ORF Size:	963 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219227).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001039091.1</a> , <a href="#">NP_001034180.1</a>
RefSeq Size:	2518 bp
RefSeq ORF:	966 bp
Locus ID:	5634
UniProt ID:	<a href="#">P11908</a>
Cytogenetics:	Xp22.2
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Pentose phosphate pathway, Purine metabolism



[View online »](#)

**MW:** 34.9 kDa

**Gene Summary:** This gene encodes a phosphoribosyl pyrophosphate synthetase that plays a central role in the synthesis of purines and pyrimidines. The encoded protein catalyzes the synthesis of 5-phosphoribosyl 1-pyrophosphate from ATP and D-ribose 5-phosphate. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]