

## Product datasheet for **RC219237L2V**

### **FBXO44 (NM\_183413) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | FBXO44 (NM_183413) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | FBXO44   |
| Synonyms:                 | FBG3; FBX6A; FBX30; Fbx44; Fbxo6a  |
| Mammalian Cell Selection: | None   |
| Vector:                   | pLenti-C-mGFP (PS100071)   |
| Tag:                      | mGFP   |
| ACCN:                     | NM_183413  |
| ORF Size:                 | 672 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC219237).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_183413.2</a> , <a href="#">NP_904320.1</a>  |
| RefSeq Size:              | 2823 bp  |
| RefSeq ORF:               | 675 bp   |
| Locus ID:                 | 93611  |
| UniProt ID:               | <a href="#">Q9H4M3</a>   |
| Cytogenetics:             | 1p36.22  |
| Protein Families:         | Druggable Genome   |
| MW:                       | 25.7 kDa   |



[View online »](#)

**Gene Summary:**

This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbxs class. It is also a member of the NFB42 (neural F Box 42 kDa) family, similar to F-box only protein 2 and F-box only protein 6. Several alternatively spliced transcript variants encoding two distinct isoforms have been found for this gene. [provided by RefSeq, Feb 2015]