

Product datasheet for **RC219270L1V**

ATG16L2 (NM_033388) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ATG16L2 (NM_033388) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ATG16L2
Synonyms:	ATG16B; WDR80
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_033388
ORF Size:	1857 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219270).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_033388.1
RefSeq Size:	2160 bp
RefSeq ORF:	1860 bp
Locus ID:	89849
UniProt ID:	Q8NAA4
Cytogenetics:	11q13.4
MW:	68.8 kDa
Gene Summary:	May play a role in autophagy.[UniProtKB/Swiss-Prot Function]


[View online »](#)