

Product datasheet for RC219320

OR56B1 (NM_001005180) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | OR56B1 (NM_001005180) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | OR56B1 |
| Synonyms: | OR11-65; OR56B1P |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| ORF Nucleotide Sequence: | >RC219320 representing NM_001005180 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATCATATGTCTGCATCTCTCAAAATCTCCAATAGTTCCAAATCCAGGTCTCTGAGTTCATCCTGC
TGGGATCCCAGGATTACAGCTGGCAACTGGCTATCTCTGCCCTGGCACTACTGTATCTCTCAGC
ACTTGCTGCAAACACCCTCATCCTCATCATCTGGCAGAACCTTCTTTACAGCAGCCCATGTATATT
TTCCTTGGCATCCTCTGTATGGTAGACATGGGTCTGGCCACTACTATCATCCCTAAGATCCTGGCCATCT
TCTGGTTTGTATGCCAAGTTATTAGCCTCCCTGAGTGCTTTGCTCAGATTTATGCCATTCACCTTCTTGT
GGGCATGGAGTCTGGTATCCTCCTCTGCATGGCTTTTGTATAGATATGTGGCTATTTGTCACCTCTTCGC
TATCCATCAATTGTCACCAAGTTCCTTAATCTTAAAAGCTACCCTGTTCATGGTGCTGAGAAATGGCTTAT
TTGCACTCCAGTGCCTGTGCTTGCAGCACAGCGTGATTATTGCTCCAAGAATGAAATTGAACACTGCCT
GTGCTCTAACCTTGGGGTCAAGCCTGGCTTGTGATGACAGGAGGCCAAACAGCATTGGCCAGTTGGTT
CTGGCATGGCTTGGAAATGGGGAGTGATCTAAGTCTTATTACTGTATATTTTATTGATTCTGTACTCTG
TACTTAGACTGAAGTCAAGCTGAAGCTGCAGCAAGGCCCTGAGCACTTGTAGTTCACATCTACCCCTCAT
CCTTTTCTTTTACACTATTGTTGTAGTGATTTTCAGTGACTCATCTGACAGAGATGAAGGCTACTTTGATT
CCAGTTCTACTTAATGTGTTGCACAACATCATCCCCCTTCCCTCAACCCTACAGTTTATGCACCTCAGA
CCAAAGAACTTAGGGCAGCCTTCCAAAAGGTGCTGTTTGCCTTACAAAAGAAATAAGATCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC219320 representing NM_001005180
Red=Cloning site Green=Tags(s)

MNHMSASLKISNSSKFQVSEFILLGFPGIHSWQHWSLPLALLYLSALAANTLILIIWQNPSLQQPMYI
 FLGILCMVDMGLATTIIPKILAIWFWDKVISLPECFAQIYAIHFFVGMESGILLCMAFDRYVAICHPLR
 YPSIVTSSLILKATLFMVLNRGLFVTPVPVLAQRDYCSKNEIEHCLCSNLGVTSLACDDRRPNISICQLV
 LAWLGMGSDLSLIILSYILILYSVLRRLNSAEAAKALSTCSSHLTLILFFYTIIVVVISVTHLTEMKATLI
 PVLLNLVHNIIPPSLNPTVYALQTKELRAAFQKVLFFALTKEIRS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6101_e04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001005180

ORF Size: 972 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

| | |
|-------------------------------|--|
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001005180.3 |
| RefSeq Size: | 975 bp |
| RefSeq ORF: | 975 bp |
| Locus ID: | 387748 |
| UniProt ID: | Q8NGI3 |
| Cytogenetics: | 11p15.4 |
| Protein Families: | Druggable Genome, Transmembrane |
| Protein Pathways: | Olfactory transduction |
| MW: | 35.9 kDa |
| Gene Summary: | Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008] |

Product images:



Circular map for RC219320