

Product datasheet for **RC219327L1V**

FOXJ3 (NM_014947) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FOXJ3 (NM_014947) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FOXJ3
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_014947
ORF Size:	1866 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219327).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014947.3
RefSeq Size:	5246 bp
RefSeq ORF:	1869 bp
Locus ID:	22887
UniProt ID:	Q9UPW0
Cytogenetics:	1p34.2
Domains:	FH
MW:	68.8 kDa



[View online »](#)

Gene Summary:

Transcriptional activator of MEF2C involved in the regulation of adult muscle fiber type identity and skeletal muscle regeneration.[UniProtKB/Swiss-Prot Function]