

Product datasheet for **RC219332A1V**

Human RNF39 (NM_025236) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human RNF39 (NM_025236) AAV Particle
Tag:	Myc-DDK
Symbol:	RNF39
Synonyms:	FAP216; HZF; HZFW; LIRF
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC219332 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTGGTGGAGAGATTTGACGAGACTGAGACTCTGGTTGAAGAGAGAGGCAATCCCAGGAGAGGGCGCGA
AAGCGGCAAAAGTTAATGCGGGAGTCGGAGAGAAGGGCATCTACACAGCAAGCAGCAGGGCGGCCCGCC
ATCTGCGCGCTCGAAGCGGTACGGTGGTCGCGGAAGGGCGGCGTCCAGATCCTGGCTTTCCATGGAT
GCGCCCGAGCTGGGCGCGGGCTGGTGGAGCGTCTGGAGCAGCTGGCGACGTGTCTCTGTGCGGGGCT
CCTTCGAGGACCGGTGCTTCTGGCGTGCAGCACAGCTTCTGCCGCGGTGTCTGGCCCGCGCTGGGG
GACTCCGCGCGGCGACCGGACCGAGGCTTCCCCACCGCTGTCCCTGCTGCGGCGTGGCGTGTCCCGC
CGCAGCCTGAGGTCTAATGTGCGGTGGCGGTGGAGGTGCGAATCAGCCGCGAGCTGCGAGAGAAGCTGG
CTGAGCCTGGGCGCGTGCAGGAGACGCCGAGGGGGCGCATCCCCACCATGGGCTGCCTGGACCTGCC
CGGAGAGGATATGAGGAAGACATGGAGACGATTTGAAGTCCCAACATCCAAGTCATCTAATTCAGAGGAT
GATCTCCCTGAAGATTATCCAGTGGTCAAAAAATGCTTCATAGACTGACAGCCGACCTGACCCTGGACC
CTGGGACCGCACACCGCCGCTGCTCATCTCCGCCACCGCCGAGCGTACAACCTGGCCCCACAGGGAC
GCCCGCGCCCCCTGACGGCCCCAAGCGCTTCGATCAGCTCCAGCTGTGCTGGGTGCGCAGGGCTTCGGG
GCCGGCGCCACTGCTGGGAGGTGGAGACTGCGGACGCGCCTCTGCAGAGACTCTTCTGGGGAGGATG
CGGACGACGAGGAGAGCACTATGCAAGTGGGCGCGGCCGGAATCAGTGCAACGCAAGGCTGCGTAAG
GCTGTGCCCTGCGGGGGCGGTGTGGGCCGTGGAGGGCCGCGCGGCCGCTGTGGGCCCTCACGGCACCC
GAACCCACCCTGCTGGGCGGTGTTGAGCCCCGCGCGCGCATTGCGTGGACCTGGACTGGGAGCGGG
GCCGCTGGCCTTCTACGACGGCGCTCACTCGACCTGCTTTACGCCTTCAGGCGCCTGGCCCCCTGGG
GGAGCGCATCTCCCGCTGTTCTGCACCTGCGACCCTCGTGCTCCGCTCCGATTGTACCAGCGGAAAGC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC219332 protein sequence Red=Cloning site Green=Tags(s) MWWRDLTRLRLWLKREAIPEGGRKAAKVNAGVGEGKIYTASSRGGPPSARSKAVTVVAEGAASRSWLSMD APELGPGLVERLEQLATCPLCGGSFEDPVLLACEHSFCRACLARRWGTPPATGTEASPTACPCGLPCPR RSLRSNVRLAVEVRISRELREKLAEPGARAGRRRGRIPTMGCLDLPGEDMRKTWRRFEVPTSKSSNSD DLPEDYPVVKMLHRLTADLTLDPGTAHRLLISADRRSVQLAPPGTPAPPDGPKRFDQLPAVLGAQGFG AGRHCWEVETADAASCRDSSGEDADDEESHYAVGAAGESVQRKGCVRLCPAGAVWAVEGRGGRLWALTAP EPTLLGGVEPPRRIRVDLDWERGRVAFYDGRSLDLLYAFQAPGPLGERIFPLFCTCDPRAPLRIVPAES TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_025236
ORF Size:	1260 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_025236.2
RefSeq Size:	2170 bp
RefSeq ORF:	1059 bp
Locus ID:	80352
UniProt ID:	Q9H2S5
Cytogenetics:	6p22.1
MW:	45.5 kDa