

Product datasheet for **RC219332L3V**

RNF39 (NM_025236) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RNF39 (NM_025236) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RNF39
Synonyms:	FAP216; HZF; HZFW; LIRF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025236
ORF Size:	1260 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219332).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025236.2
RefSeq Size:	2170 bp
RefSeq ORF:	1059 bp
Locus ID:	80352
UniProt ID:	Q9H2S5
Cytogenetics:	6p22.1
Protein Families:	Druggable Genome
MW:	45.5 kDa



[View online »](#)

Gene Summary:

This gene lies within the major histocompatibility complex class I region on chromosome 6. Studies of a similar rat protein suggest that this gene encodes a protein that plays a role in an early phase of synaptic plasticity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]