

## Product datasheet for **RC219410L4V**

### **U1SNRNPBP (SNRNP35) (NM\_180703) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	U1SNRNPBP (SNRNP35) (NM_180703) Human Tagged ORF Clone Lentiviral Particle
Symbol:	U1SNRNPBP
Synonyms:	HM-1; MGC138160; small nuclear ribonucleoprotein 35kDa (U11/U12); U1 snRNP binding protein homolog; U1-snRNP binding protein homolog; U1SNRNPBP; U11/U12 snRNP 35K
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_180703
ORF Size:	753 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219410).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_180703.1</a> , <a href="#">NP_851034.1</a>
RefSeq Size:	1523 bp
RefSeq ORF:	755 bp
Locus ID:	11066
Cytogenetics:	12q24.31
MW:	30 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a homolog of the U1-snRNP binding protein. The N-terminal half contains a RNA recognition motif and the C-terminal half is rich in Arg/Asp and Arg/Glu dipeptides, which is a characteristic of a variety of splicing factors. This protein is a component of the U11/U12 small nuclear ribonucleoproteins (snRNP) that form part of the U12-type spliceosome. Alternative splicing results in multiple transcript variants encoding two distinct isoforms and representing a non-protein coding variant. [provided by RefSeq, Aug 2013]