

Product datasheet for **RC219423L3V**

AP4 (REPIN1) (NM_001099695) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AP4 (REPIN1) (NM_001099695) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AP4
Synonyms:	AP4; RIP60; Zfp464; ZNF464
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001099695
ORF Size:	1872 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219423).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001099695.1
RefSeq Size:	3185 bp
RefSeq ORF:	1875 bp
Locus ID:	29803
UniProt ID:	Q9BWE0
Cytogenetics:	7q36.1
MW:	69.7 kDa



[View online »](#)

Gene Summary:

Sequence-specific double-stranded DNA-binding protein required for initiation of chromosomal DNA replication. Binds on 5'-ATT-3' reiterated sequences downstream of the origin of bidirectional replication (OBR) and a second, homologous ATT sequence of opposite orientation situated within the OBR zone. Facilitates DNA bending.[UniProtKB/Swiss-Prot Function]