

Product datasheet for **RC219434L3V**

FUT8 (NM_004480) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FUT8 (NM_004480) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FUT8
Synonyms:	CDGF; CDGF1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004480
ORF Size:	1236 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219434).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004480.4 , NP_004471.4
RefSeq Size:	4749 bp
RefSeq ORF:	1239 bp
Locus ID:	2530
UniProt ID:	Q9BYC5
Cytogenetics:	14q23.3
Domains:	SH3
Protein Families:	Transmembrane



[View online »](#)

Protein Pathways: Keratan sulfate biosynthesis, Metabolic pathways, N-Glycan biosynthesis

MW: 47.3 kDa

Gene Summary: This gene encodes an enzyme belonging to the family of fucosyltransferases. The product of this gene catalyzes the transfer of fucose from GDP-fucose to N-linked type complex glycopeptides. This enzyme is distinct from other fucosyltransferases which catalyze alpha1-2, alpha1-3, and alpha1-4 fucose addition. The expression of this gene may contribute to the malignancy of cancer cells and to their invasive and metastatic capabilities. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2011]