

Product datasheet for **RC219463L3V**

RHBDD2 (NM_001040457) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | RHBDD2 (NM_001040457) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | RHBDD2 |
| Synonyms: | NPD007; RHBDL7 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001040457 |
| ORF Size: | 669 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC219463). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001040457.1 , NP_001035547.1 |
| RefSeq Size: | 1878 bp |
| RefSeq ORF: | 672 bp |
| Locus ID: | 57414 |
| UniProt ID: | Q6NTF9 |
| Cytogenetics: | 7q11.23 |
| Protein Families: | Transmembrane |
| MW: | 23.7 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene is a member of the rhomboid family of membrane-bound proteases and is overexpressed in some breast cancers. Members of this family are involved in intramembrane proteolysis. In mouse, the orthologous protein associates with the Golgi body. [provided by RefSeq, Sep 2016]