

## Product datasheet for **RC219468L3V**

### SF4 (SUGP1) (NM\_172231) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SF4 (SUGP1) (NM_172231) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SF4
Synonyms:	F23858; RBP; SF4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172231
ORF Size:	1935 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219468).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172231.2</a>
RefSeq Size:	2101 bp
RefSeq ORF:	1938 bp
Locus ID:	57794
UniProt ID:	<a href="#">Q8IWZ8</a>
Cytogenetics:	19p13.11
MW:	72.3 kDa
Gene Summary:	SF4 is a member of the SURP family of splicing factors.[supplied by OMIM, Sep 2003]


[View online »](#)