




<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_153046.2</a> , <a href="#">NP_694591.2</a>
<b>RefSeq Size:</b>	4782 bp
<b>RefSeq ORF:</b>	4149 bp
<b>Locus ID:</b>	122402
<b>UniProt ID:</b>	<a href="#">Q8NDG6</a>
<b>Cytogenetics:</b>	14q32.33
<b>MW:</b>	155.5 kDa
<b>Gene Summary:</b>	ATP-binding RNA helicase required during spermatogenesis (PubMed:28536242). Required to repress transposable elements and prevent their mobilization, which is essential for the germline integrity. Acts via the piRNA metabolic process, which mediates the repression of transposable elements during meiosis by forming complexes composed of piRNAs and Piwi proteins and governs the methylation and subsequent repression of transposons. Acts downstream of piRNA biogenesis: exclusively required for transposon silencing in the nucleus, suggesting that it acts as a nuclear effector in the nucleus together with PIWIL4. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC219472L3