

Product datasheet for **RC219475L4V**

TRAF3IP3 (NM_025228) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TRAF3IP3 (NM_025228) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRAF3IP3
Synonyms:	T3JAM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025228
ORF Size:	1653 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219475).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025228.1
RefSeq Size:	2058 bp
RefSeq ORF:	1656 bp
Locus ID:	80342
UniProt ID:	Q9Y228
Cytogenetics:	1q32.2
Protein Families:	Druggable Genome, Transmembrane
MW:	63.6 kDa



[View online »](#)

Gene Summary:

The gene encodes a protein that mediates cell growth by modulating the c-Jun N-terminal kinase signal transduction pathway. The encoded protein may also interact with a large multi-protein assembly containing the phosphatase 2A catalytic subunit. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013]