

Product datasheet for **RC219478A1V**

Human NCKAP5L (NM_001037806) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human NCKAP5L (NM_001037806) AAV Particle
Tag: Myc-DDK
Symbol: NCKAP5L
Synonyms: Cep169; FP1193; KIAA1602
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC219478 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGACCAGCCAGCTGGGGTCTGGAACCCAAGGCCAGGAGAGGGTGATGATGGCAGCATGGAGCCAG
GCACCTGCCAGGAGCTTCTGCACCGACTGCGGGAGCTGGAGGCAGAGAACTCGGCACTTGCCAGGCCAA
CGAAAACAGCGGGAGACTTATGAGCGCTGTCTGGACGAGGTCTGTGGGTCTGTAGTGGGACTGGGGGA
TGTGGCTCATCTGCTCCTGGCAGAAGCTGGGGTCAGCTGATGGCTCTGCCTCGGGGCTTTCTGTCCCCAG
GTTGCCAACCATGTGGTACAGGCGTTGCTGAACCAGAAGGTGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC219478 protein sequence
Red=Cloning site Green=Tags(s)
MDQPAGGPGNPRPGEEDGSMEPGTCQELLHRLRELEAENSALAQANENQRETYERCLDEVCGSVVGLGG
CGSSAPGRSWGQLMALPRGFLSPGCQPCGTGVAEPEGE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001037806
ORF Size: 324 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001037806.1</u>
RefSeq Size:	4918 bp
RefSeq ORF:	4005 bp
Locus ID:	57701
UniProt ID:	<u>Q9HCH0</u>
Cytogenetics:	12q13.12
MW:	11.2 kDa