

Product datasheet for **RC219522**

PEX5 related protein (PEX5L) (NM_016559) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PEX5 related protein (PEX5L) (NM_016559) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PEX5 related protein
Synonyms:	PEX5R; PEX5RP; PXR2; PXR2B; TRIP8b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC219522 representing NM_016559
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTACCAGGGACACATGCAGAAAAGTAAAGAACAAGGATATGGAAACTAAGCAGTGATGAAGACCTCG
 AAATAATTGTTGATCAAAAGCAGGGAAAAGGCTCTAGGGCGGCAGATAAGGCTGTTGCCATGGTGATGAA
 GGAGATACCGAGGGAGGAGTCTGCTGAAGAAAAGCCCTCCTTACTATGACATCACAGCTGGTGAATGAG
 CAACAAGAAAGCAGACCCCTCTGAGTCCCTCCATCGATGACTTTCTCTGTGAAACCAAAATCGGAAGCAA
 TAGCAAGGCCAGTAACATCCAATACAGCTGTATTGACCACTGGCTTAGATCTCCTCGACCTGAGTGAACC
 AGTCTCTCAAACCCAAACCAAGCCAAGAAGTCAGAGCCCTCATCAAAACCTCATCCCTCAAGAAAAAG
 GCCGATGGATCTGACCTCATCAGCACGGATGCTGAGCAGAGAGGCCAGCCTCTCAGAGTCCCGGAGACTT
 CATCCTTAGATCTAGACATTCAACACAACCTGGAAAAATGGGACGATGTTAAGTTTCATGGAGATCGAAA
 TACCAAGGGACATCCAATGGCAGAGAGAAAATCATCCTCATCTAGAACTGGATCAAAAGAGCTCTTATGG
 TCCTCAGAACACAGATCTCAACCAGAACTGAGTGGTGGAAAAAGCGCCCTCAACTCTGAGTCGGCTTCAG
 AATTGGAATTAGTGGCTCCGACTCAGGCTCGACTGACCAAGAACATCGTGGGGAAGCGCATTACTTTTC
 TAGAAACCACTCCTTAGAAGAAGAGTTTGAAAGGGCAAAAGCAGCAGTGGAGTCAGATACAGAGTTTGG
 GATAAGATGCAAGCAGAAATGGGAAGAAATGGCTCGGAGGAACTGGATATCTGAGAACCAAGAAGCCGAGA
 ACCAAGTAACCATCTCGGCTAGTGAGAAGGGATATTACTTTCACACTGAAAACCCCTTCAAGGACTGGCC
 TGGAGCATTGGAAGAAGGCTTAAAAAGGCTGAAGGAAGGGGATCTGCCAGTCACCATCCTGTTTCATGGAA
 GCAGCAATTCTTCAGGACCTGGAGATGCAGAGGCATGGCAGTTCTCGGGATAACCCAGGCGGAGAATG
 AAAATGAACAAGCAGCTATTGTGCCCCTCCAGAGGTGCTTAGAATTACAGCCCAAACTTAAAAAGCTTT
 GATGGCCTTGGCTGTGAGTTATACTAACACTGGCCATCAGCAGGATGCCTGTGACGCTCTGAAGAATTGG
 ATTAAGCAAAATCCAAGTACAATACCTTGTAAGGAAGCAAGAAGGGATCTCCAGGCCTCACCCGCGGA
 TGTCTAAGTCCCCAGTTGATAGCTCTGTTCTGGAAGGGGTGAAGGAATTATCTGGAAGCTGCCACCA
 AAATGGAGATATGATCGACCCAGACCTGCAGACAGGTCTAGGGGTTCTGTTCCACCTGAGTGGAGAATTT
 AATAGAGCAATAGATGCATTTAACGCTGCCTTAAGTGTTCGGCCAGAGGACTATTCATATGGAACCGCC
 TCGGGGCGACCTTGGCGAACGGAGACCGCAGCGAGGAAGCCGTGGAGGCCTATACGCGAGCACTGGAGAT
 TCAGCCAGGATTCATCCGGTCCAGATACAACCTAGGAATAAGCTGCATCAACCTGGGCGCCTACAGAGAA
 GCGGTGAGCAATTTCTCACTGCCCTCAGTTTGCAAAGAAAGAGCAGGAATCAGCAGCAAGTTCTCATC
 CTGCAATCTCTGGAATATCTGGGCTGCCCTCAGAATTGCGCTCTCTCTGATGGACCAACCAAGAACTCTT
 CCAGGCGGCTAATCTTGGTGACCTGGATGTCCTCTTAAGAGCTTTCAACTTGGATCCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC219522 representing NM_016559
 Red=Cloning site Green=Tags(s)

MYQGHMQKSKEQGYGKLSSDEDLIIIVDQKQGKGSRAADKAVAMVMKEIPREESAEEKPLLTMTSQLVNE
 QQESRPLLSPIIDDFLCETKSEAIARPVTSNTAVLTTGLDLLDLSEPVSTQTKAKKSEPSSTKSSLKKK
 ADGSDLISTDAEQRGQPLRVPETSSLDLDIQTQLEKWDDVKFHGDRNTKGHPMAERKSSSRSGSKELLW
 SSEHRSQPELSSGKSALNSESALEL VAPTQARLTKEHRWGSALLSRNHSLEEEFERAKAAVESDTEFW
 DKMQAEWEEMARRNWIENQEAQNQVTISASEKGYFHTENPFKDWPGAFEEGLKRLKEGDLPTVILFME
 AAILQDPGDAEAWQFLGITQAENENEQAAIVALQRCLQLPNNLKALMALAVSYTNTGHQQDACDALKNW
 IKQNPKYKYLKSKKSGPLTRMSKSPVDSSVLEGVKELYLEAAHQNGDMIDPDLQTLGVL FHLSGEF
 NRAIDAFNAALTVPEDYSLWNRLGATLANGDRSEEAVEAYTRALEIQPGFIRSRYNLGISCINLGAYRE
 AVSNFLTALSLQRKSRNQVPHPAISGNIWAALRIALSLMDQPELFQAANLGDLVDVLLRAFNLDP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM 016559

ORF Size: 1878 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: NM 016559.3

RefSeq Size: 3577 bp


RefSeq ORF: 1881 bp

Locus ID:	51555
UniProt ID:	Q8IYB4
Cytogenetics:	3q26.33
Domains:	TPR
Protein Families:	Druggable Genome
MW:	69.5 kDa
Gene Summary:	Accessory subunit of hyperpolarization-activated cyclic nucleotide-gated (HCN) channels, regulating their cell-surface expression and cyclic nucleotide dependence.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC219522



Coomassie blue staining of purified PEX5L protein (Cat# [TP319522]). The protein was produced from HEK293T cells transfected with PEX5L cDNA clone (Cat# RC219522) using MegaTran 2.0 (Cat# [TT210002]).