

## Product datasheet for **RC219526L3V**

### HPS1 (NM\_000195) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HPS1 (NM_000195) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HPS1
Synonyms:	BLOC3S1; HPS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000195
ORF Size:	2100 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219526).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000195.2</a>
RefSeq Size:	3714 bp
RefSeq ORF:	2103 bp
Locus ID:	3257
UniProt ID:	<a href="#">Q92902</a>
Cytogenetics:	10q24.2
MW:	79.1 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a protein that may play a role in organelle biogenesis associated with melanosomes, platelet dense granules, and lysosomes. The encoded protein is a component of three different protein complexes termed biogenesis of lysosome-related organelles complex (BLOC)-3, BLOC4, and BLOC5. Mutations in this gene are associated with Hermansky-Pudlak syndrome type 1. Alternative splicing results in multiple transcript variants. A pseudogene related to this gene is located on chromosome 22. [provided by RefSeq, Aug 2015]