

Product datasheet for **RC219531L4V**

ZNF117 (NM_015852) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF117 (NM_015852) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF117
Synonyms:	H-plk; HPPF9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015852
ORF Size:	1449 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219531).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015852.3
RefSeq Size:	6404 bp
RefSeq ORF:	1452 bp
Locus ID:	51351
UniProt ID:	Q03924
Cytogenetics:	7q11.21
Domains:	zf-C2H2
MW:	56.8 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein containing multiple C2H2-type zinc finger motifs. Readthrough transcription occurs between this gene and the upstream endogenous retrovirus group 3 member 1 (ERV3-1) locus, and may result in additional transcript variants encoding the zinc finger protein. [provided by RefSeq, Jan 2017]