

Product datasheet for **RC219546L3V**

ZNRF3 (NM_032173) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNRF3 (NM_032173) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNRF3
Synonyms:	BK747E2.3; RNF203
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032173
ORF Size:	2508 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219546).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032173.2 , NP_115549.2
RefSeq Size:	6541 bp
RefSeq ORF:	2511 bp
Locus ID:	84133
UniProt ID:	Q9ULT6
Cytogenetics:	22q12.1
Protein Families:	Druggable Genome, Transmembrane
MW:	89.7 kDa



[View online »](#)

Gene Summary:

E3 ubiquitin-protein ligase that acts as a negative regulator of the Wnt signaling pathway by mediating the ubiquitination and subsequent degradation of Wnt receptor complex components Frizzled and LRP6. Acts on both canonical and non-canonical Wnt signaling pathway. Acts as a tumor suppressor in the intestinal stem cell zone by inhibiting the Wnt signaling pathway, thereby restricting the size of the intestinal stem cell zone (PubMed:22575959). Along with RSPO2 and RNF43, constitutes a master switch that governs limb specification (By similarity).[UniProtKB/Swiss-Prot Function]