

Product datasheet for **RC219564L3V**

NME2 (NM_002512) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NME2 (NM_002512) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NME2
Synonyms:	NDKB; NDPK-B; NDPKB; NM23-H2; NM23B; PUF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002512
ORF Size:	456 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219564).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002512.2
RefSeq Size:	822 bp
RefSeq ORF:	459 bp
Locus ID:	4831
UniProt ID:	P22392
Cytogenetics:	17q21.33
Domains:	NDK
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

Protein Pathways: Metabolic pathways, Purine metabolism, Pyrimidine metabolism

MW: 17.1 kDa

Gene Summary: Nucleoside diphosphate kinase (NDK) exists as a hexamer composed of 'A' (encoded by NME1) and 'B' (encoded by this gene) isoforms. Multiple alternatively spliced transcript variants have been found for this gene. Read-through transcription from the neighboring upstream gene (NME1) generates naturally-occurring transcripts (NME1-NME2) that encode a fusion protein comprised of sequence sharing identity with each individual gene product. [provided by RefSeq, Nov 2010]