

Product datasheet for **RC219569A1V**

Human CD32B (FCGR2B) (NM_001002273) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CD32B (FCGR2B) (NM_001002273) AAV Particle
Tag:	Myc-DDK
Symbol:	CD32B
Synonyms:	CD32; CD32B; FCG2; FCGR2; FCGR2C; FcRII-c; IGFR2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC219569 representing NM_001002273 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGAATCCTGTCATTCTTACCTGTCCTTGCCACTGAGAGTGACTGGGCTGACTGCAAGTCCCCCAGC
CTTGGGGTCATATGCTTCTGTGGACAGCTGTGCTATTCCTGGCTCCTGTTGCTGGGACACCTGCTCCCC
AAAGGCTGTGCTGAACTCGAGCCCCAGTGGATCAACGTGCTCCAGGAGGACTCTGTGACTCTGACATGC
CGGGGGACTCACAGCCCTGAGAGCGACTCCATTCACTGGTTCACAATGGGAATCTCATTCCACCCACA
CGCAGCCAGCTACAGTTCAAGGCCAACAAATGACAGCGGGGAGTACACGTGCCAGACTGGCCAGAC
CAGCCTCAGCGACCCTGTGCATCTGACTGTGCTTTCTGAGTGGCTGGTGTCCAGACCCCTCACCTGGAG
TTCCAGGAGGGAGAAACCATCGTGCTGAGGTGCCACAGCTGGAAGGACAAGCCTCTGGTCAAGGTCACAT
TCTTCCAGAATGGAAAATCCAAGAAATTTCCCGTTCCGATCCCAACTTCTCCATCCCACAAGCAAACCA
CAGTCACAGTGGTGATTACCACTGCACAGGAACATAGGCTACACGCTGTACTCATCCAAGCCTGTGACC
ATCACTGTCCAAGCTCCAGCTTTCACCGATGGGGATCATTGTGGCTGTGGTCACTGGGATTGCTGTAG
CGGCCATTGTTGCTGCTGTAGTGGCCTTGATCTACTGCAGGAAAAAGCGGATTTACGCCAATCCCACTAA
TCCTGATGAGGCTGACAAAGTTGGGGCTGAGAACACAATCACCTATTCATTCTCATGCACCCGGATGCT
CTGGAAGAGCCTGATGACCAGAACCGTATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC219569 representing NM_001002273 Red=Cloning site Green=Tags(s) MGILSFLPVLATESDWADCKSPQPWGHMLLWTAVLFLAPVAGTPAPPKAVLKLEPQWINVLQEDSVTLTC RGTHSPESDSIQWFHNGNLIPTHTQPSYRFKANNNDSGEYTCQTGQTSLSDPVHLTVLSEWLVLQTPHLE FQEGETIVLRCHSWKDKPLVKVTFQNGKSKKFSRSDPNFSIPQANHSHSGDYHCTGNIGYTLYSSKPVT ITVQAPSSSPMGIIVAVVTGIAVAAIVAAVVALIYCRKKRISANPTNPDEADKVGAEITYSLLMHPDA LEEDDQNR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001002273
ORF Size:	870 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001002273.1</u>
RefSeq Size:	1573 bp
RefSeq ORF:	873 bp
Locus ID:	2213
UniProt ID:	<u>P31994</u>
Cytogenetics:	1q23.3
MW:	31.7 kDa