

## Product datasheet for **RC219571L2V**

### ZNF14 (NM\_021030) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF14 (NM_021030) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF14
Synonyms:	GIOT-4; KOX6
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_021030
ORF Size:	1926 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219571).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021030.2</a>
RefSeq Size:	2967 bp
RefSeq ORF:	1929 bp
Locus ID:	7561
UniProt ID:	<a href="#">P17017</a>
Cytogenetics:	19p13.11
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 75.2 kDa

**Gene Summary:** The protein encoded by this gene contains a zinc finger and a Kruppel-associated box (KRAB) domain. KRAB domain is known to be involved in the transcriptional repression of a number of zinc finger proteins. [provided by RefSeq, Jul 2008]