

Product datasheet for **RC219636L4V**

TMEM180 (MFSD13A) (NM_024789) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM180 (MFSD13A) (NM_024789) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM180
Synonyms:	bA18114.8; C10orf77; TMEM180
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024789
ORF Size:	1551 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219636).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024789.2 , NP_079065.2
RefSeq Size:	2846 bp
RefSeq ORF:	1554 bp
Locus ID:	79847
UniProt ID:	Q14CX5
Cytogenetics:	10q24.32
Protein Families:	Transmembrane
MW:	57.2 kDa



[View online »](#)