

Product datasheet for RC219639A1V

Human PATE (PATE1) (NM_138294) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PATE (PATE1) (NM_138294) AAV Particle
Tag: Myc-DDK
Symbol: PATE
Synonyms: PATE
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC219639 representing NM_138294
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGACAAGTCCCTCTTGCTGGAACCCCCATCCTGCTGCTGCTTTAGGGCATTATCTGGATCACTTT
 CAATGAGAAATGATGCAGTCAATGAAATAGTTGCTGTGAAAAACAATTTCTGTGATAGAAATTGTTCA
 GTGTAGGATGTGCCACCTCCAGTCCAGGAGAAAAGTGCTCCAGAGGAAGAGGAATATGCACAGCAACA
 ACAGAAGAGGCCTGCATGGTTGGAAGGATGTTCAAAGGGATGGTAATCCCTGGTTAACCTTCATGGGCT
 GCCTAAAGAACTGTGCTGATGTGAAAGGCATAAGGTGGAGTGTCTATTTGGTGAACCTTCAGGTGCTGCAG
 GAGCCATGACCTGTGCAATGAAGACCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC219639 representing NM_138294
 Red=Cloning site Green=Tags(s)

MDKSLLELPILLCCFRALSGSLSMRNDVNEIVAVKNNFPVIEIVQCRMCHLQFPGEKCSRGRGICTAT
 TEEACMVGMRMFKRDGNPWLTFMGCLKNCADVKGIRWSVYLVNFRCCRS HDLCNEDL

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_138294


[View online »](#)

ORF Size:	378 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_138294.2</u>
RefSeq Size:	1576 bp
RefSeq ORF:	381 bp
Locus ID:	160065
UniProt ID:	<u>Q8WXA2</u>
Cytogenetics:	11q24.2
MW:	14.1 kDa