

Product datasheet for **RC219650L3V**

FAM188B (MINDY4) (NM_032222) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM188B (MINDY4) (NM_032222) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MINDY4
Synonyms:	AQP-1; AQP1; C7orf67; CHIP28; FAM188B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032222
ORF Size:	2271 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219650).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032222.2
RefSeq Size:	2750 bp
RefSeq ORF:	2274 bp
Locus ID:	84182
UniProt ID:	Q4G0A6
Cytogenetics:	7p14.3
MW:	84.4 kDa


[View online »](#)

Gene Summary: Probable hydrolase that can remove 'Lys-48'-linked conjugated ubiquitin from proteins.
[UniProtKB/Swiss-Prot Function]