

Product datasheet for **RC219659**

Heparanase 1 (HPSE) (NM_001098540) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Heparanase 1 (HPSE) (NM_001098540) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Heparanase 1
Synonyms:	HPA; HPA1; HPR1; HPSE1; HSE1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC219659 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTGCTGCGCTCGAAGCCTGCGCTGCCGCCCGCTGATGCTGCTGCTCCTGGGGCCGCTGGGTCCCC
 TCTCCCTGGCGCCTGCCCGACTGCGCAAGCACAGGACGTCGTGGACCTGGACTTCTCACCCAGGA
 GCCGCTGCACCTGGTGGAGCCCTCGTTCCTGTCGGTACACATTGACGCCAACCTGGCCACGACCCGCGG
 TTCCTCATCCTCTGGGTTCTCCAAAGCTTCGTACCTGGCCAGAGGCTTGTCTCTGCGTACCTGAGGT
 TTGGTGGACCAAGACAGACTTCTAATTTTCGATCCCAAGAAGGAATCAACCTTTGAAGAGAGAAGTTA
 CTGGCAATCTCAAGTCAACCAGGATATTGCAAAATATGGATCCATCCCTCCTGATGTGGAGGAGAAGTTA
 CGGTTGGAATGGCCCTACCAGGAGCAATTGCTACTCCGAGAACACTACCAGAAAAAGTTCAAGAACAGCA
 CCTACTCAAGAAGCTCTGTAGATGTGCTATACACTTTTGCAAAGCTGCTCAGGACTGGACTTGATCTTTGG
 CCTAAATGCGTTATTAAGAACAGCAGATTTGCAGTGGAAACAGTTCTAATGCTCAGTTGCTCTGGACTAC
 TGCTCTTCCAAGGGGTATAACATTTCTGGGAAGTGGCAATGAACCTAACAGTTTCTTAAAGAAGGCTG
 ATATTTTCATCAATGGTTCGAGTTAGGAGAAGATTTTATTCAATTGCATAAACTCTAAGAAAAGTCCAC
 CTTCAAAAATGCAAACTCTATGGTCTGATGTTGGTCAGCCTCGAAGAAAGACGGCTAAGATGCTGAAG
 AGCTTCTGAAGGCTGGTGGAGAAGTATTGATTCAGTTACATGGCATCACTACTATTTGAATGGACGGA
 CTGCTACCAGGGAAGATTTTCTAAACCTGATGTATTGGACATTTTTATTTATCTGTGCAAAAAGTTTT
 CCAGGTGGTTGAGAGCACCAGGCTGGCAAGAAGTCTGGTTAGGAGAAACAAGCTCTGCATATGGAGGC
 GGAGGCCCTTGTATCCGACACCTTTGAGCTGGCTTTATGTGGCTGGATAAATGGGCTGTCAGCCC
 GAATGGGAATAGAAGTGGTATGAGGCAAGTATTCTTTGGAGCAGAACTACCATTAGTGGATGAAAA
 CTTGATCCTTTACCTGATTATTGGCTATCTCTTCTGTTCAAGAAATTGGTGGCACCAAGGTGTTAATG
 GCAAGCGTGCAAGGTTCAAGAGAAGGAAGCTTCGAGTATACCTTCATTGCACAAACACTGACAATCCAA
 GGTATAAGAAGGAGATTTAACTCTGTATGCCATAAACCTCCATAATGTCACCAAGTACTTGCAGTTACC
 CTATCCTTTTTCTAACAAGCAAGTGGATAAATACCTTCTAAGACCTTTGGGACCTCATGGATTACTTTCC
 AAATCTGTCCAACCTCAATGGTCTAACTCTAAAGATGGTGGATGATCAAACCTTGCACCTTTAATGGAAA
 AACCTCTCGGCCAGGAAGTTCAGTGGCTTGGCAGCTTCTCATATAGTTTTTTTGTGATAAGAAATGC
 CAAAGTTGCTGCTGCATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC219659 protein sequence
 Red=Cloning site Green=Tags(s)

MLLRSPALPPPLMLLLLGPLGPLSPGALPRPAQAQDVVDLDFFTQEPLHLVSPSFLSVTIDANLATDPR
 FLILLGSPKLRTLARGLSPAYLRFGGTKTDFLIFDPKKESTFEERSYWQSQVNQDICKYGSIPPDVEEKL
 RLEWPYQEQLLLREHYQKFKNSTYSRSSVDVLYTFANCSGLDLIFGLNALLRTADLQWNSSNAQLLLDY
 CSSKGYNISWELGNEPNSFLKKADIFINGSQLGEDFIQLHKLKSTFKNAKLYGPDVGGPQRRTAKMLK
 SFLKAGGEVIDSVTWHHYLNGRTATREDFLNPDVLDIFISSVQKVFQVVESTRPKQVWLGETSSAYGG
 GAPLLSDTFAAGFMWLDKLGLSARMGIEVVMRQVFFGAGNYHLVDENFDPLPDYWL SLLFKKLVGTVKVL
 ASVQGSRRKLRVYLHCTNTDNPRYKEDLTLYAINLHNVTKYLRLPYPF SNKQVDKYLRLPLGPHGLLS
 KSVQLNGLTLKMVDDQTLPLMEKPLRPGSSLGLPAFSYSFFVIRNAKVAACI

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6149_e12.zip

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001098540

ORF Size: 1629 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001098540.3](#)

RefSeq Size: 4630 bp

RefSeq ORF: 1632 bp

Locus ID: 10855

UniProt ID: [Q9Y251](#)

Cytogenetics: 4q21.23


Protein Families: Secreted Protein

Protein Pathways: Glycosaminoglycan degradation, Metabolic pathways


MW: 61.2 kDa

Gene Summary: Heparan sulfate proteoglycans are major components of the basement membrane and extracellular matrix. The protein encoded by this gene is an enzyme that cleaves heparan sulfate proteoglycans to permit cell movement through remodeling of the extracellular matrix. In addition, this cleavage can release bioactive molecules from the extracellular matrix. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]

Product images:



Circular map for RC219659



Western blot validation of overexpression lysate (Cat# [LY420636]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC219659 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).