

Product datasheet for **RC219670L3V**

ZNF324 (NM_014347) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF324 (NM_014347) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF324
Synonyms:	ZF5128; ZNF324A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014347
ORF Size:	1659 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219670).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014347.1 , NP_055162.1
RefSeq Size:	3000 bp
RefSeq ORF:	1662 bp
Locus ID:	25799
UniProt ID:	O75467
Cytogenetics:	19q13.43
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors


[View online »](#)

MW: 60.9 kDa

Gene Summary: May be involved in transcriptional regulation. May be involved in regulation of cell proliferation.[UniProtKB/Swiss-Prot Function]