

## Product datasheet for RC219681

### CYLC2 (NM\_001340) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** CYLC2 (NM\_001340) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** CYLC2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC219681 representing NM\_001340  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCGTTATGCTCCATGTATTCTGCAACTTTTCATCATACATTAGGTTTTGGCGATTTAGCCATAATTT  
 GGCATGAATTCAGGTCTTTTAAGTTTTATTCCATTGTAAGAATAAACAAGTTTGTTTCATTCATGTCTCCG  
 AAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC219681 representing NM\_001340  
 Red=Cloning site Green=Tags(s)  
 MALCSMYSATFHHTLGFGLAIIWHEFRSFKFYIVRINKFVHSLRKKKKKKKKK

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI


[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001340

**ORF Size:** 168 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<u>NM_001340.1, NP_001331.1</u>
<b>RefSeq Size:</b>	2119 bp
<b>RefSeq ORF:</b>	1047 bp
<b>Locus ID:</b>	1539
<b>UniProt ID:</b>	<u>Q14093</u>
<b>Cytogenetics:</b>	9q31.1
<b>MW:</b>	6.74 kDa
<b>Gene Summary:</b>	Cylicin II (CYLC2) is specifically expressed in testis and is part of the cytoskeletal calyx of mammalian sperm heads. Cylicin II may play a role in the morphogenesis of the sperm head. [provided by RefSeq, Jul 2008]

## Product images:



Circular map for RC219681