

## Product datasheet for **RC219696L2V**

### TMEM26 (NM\_178505) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM26 (NM_178505) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM26
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_178505
ORF Size:	1104 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219696).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178505.5</a> , <a href="#">NP_848600.2</a>
RefSeq Size:	5170 bp
RefSeq ORF:	1107 bp
Locus ID:	219623
UniProt ID:	<a href="#">Q6ZUK4</a>
Cytogenetics:	10q21.2
Protein Families:	Transmembrane
MW:	41.5 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein containing multiple transmembrane helices. It is a selective surface protein marker of brite/beige adipocytes, which may coexist with classical brown adipocytes in brown adipose tissue. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Dec 2015]