

Product datasheet for **RC219761L2V**

ADAR1 (ADAR) (NM_001025107) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ADAR1 (ADAR) (NM_001025107) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADAR1
Synonyms:	ADAR1; AGS6; DRADA; DSH; DSRAD; G1P1; IFI-4; IFI4; K88DSRBP; P136
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001025107
ORF Size:	2796 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219761).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001025107.3
RefSeq Size:	6532 bp
RefSeq ORF:	2796 bp
Locus ID:	103
UniProt ID:	P55265
Cytogenetics:	1q21.3
Protein Families:	Druggable Genome
Protein Pathways:	Cytosolic DNA-sensing pathway



[View online »](#)

MW: 103.6 kDa

Gene Summary: This gene encodes the enzyme responsible for RNA editing by site-specific deamination of adenosines. This enzyme destabilizes double-stranded RNA through conversion of adenosine to inosine. Mutations in this gene have been associated with dyschromatosis symmetrica hereditaria. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2010]