

Product datasheet for **RC219778A1V**

Human **OR9G1** (NM_001005213) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human OR9G1 (NM_001005213) AAV Particle
Tag:	Myc-DDK
Symbol:	OR9G1
Synonyms:	OR9G5
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC219778 representing NM_001005213 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAGAGGAGCAATCATACAGTGACTGAGTTTATACTGCTGGGCTTCACCACAGACCCAGGAATGCAGC
TGGGCTCTTCGTGGTGTTCCTGGGCGTGTACTCTCTCACTGTGGTAGGAAATAGCACCTCATCGTGT
GATCTGTAATGACTCCTGCCTCCACACACCCATGTATTTTCACTGGAAATCTGTCGTTTCTGGATCTC
TGGTATTCTTCTGTCTACACCCAAAGATCCTAGTGACCTGCATCTCTGAAGACAAAAGCATCTCCTTTG
CTGGCTGCCTGTGTGAGTCTTCTTCTCTGCAGGCTGGCCTATAGTGAGTGCTACCTGCTGGCTGCCGT
GGCTTATGACCGCTACGTGGCCATCTCCAAGCCCTGCTTTATGCCAGGCCATGTCCATAAAGCTGTGT
GCATTGCTGGTAGCAGTCTCATATTGTGGTGGCTTTATTAACCTTTCAATCATCACCAAGAAAACGTTTT
CCTTTAACTTCTGCCGTGAAAACATCATTGATGACTTTTTCTGTGATTTGCTTCCTTGGTGGAGCTGGC
CTGTGGCGAGAAGGGCGGCTATAAAATTATGATGTACTTCTGCTGGCCTCCAATGTCATCTGCCCCGCA
GTGCTCATCTGGCCTCCTACCTCTTTATCATCACCAGTGTCTTGAGGATCTCCTCCTCAAGGGCTACC
TCAAAGCCTTCTCCACATGCTCCTCCACCTGACCTCTGTCACTTTATACTATGGCTCCATTCTCTACAT
CTACGCTCTCCCAGATCTAGCTATTCTTTGATATGGACAAAATAGTTTCTACATTTTACACTGTGGTA
TTCCCCATGTTGAATCTCATGATCTACAGCCTAAGGAATAAGGATGTGAAAGAGGCTCTGAAAAAATTC
TCCCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC219778 representing NM_001005213 Red=Cloning site Green=Tags(s) MQRSNHTVTEFILLGFTTDPGMQLGLFVVFLGVYSLTVVGNSTLIVLICNDSCLHTPMYFFTGNLSFLDL WYSSVYTPKILVTCISEDKSISFAGCLCQFFFAGLAYSECYLLAAVAYDRYVAISKPLLYAQAMSIKLC ALLVAVSYCGGFINSIIITKKTFSFNFCRENIIDDFCDLLPLVELACGEKGGYKIMMYFLLASNVICPA VLILASYLFIITSVLRISSSKGYLKAFSTCSSHLTSVTLYYGSILYIYALPRSSYSFDMDKIVSTFYTVV FPMLNLMYSLRNKDVKEALKLLP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001005213
ORF Size:	915 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001005213.1, NP_001005213.1</u>
RefSeq Size:	918 bp
RefSeq ORF:	918 bp
Locus ID:	390174
UniProt ID:	<u>Q8NH87</u>
Cytogenetics:	11q12.1
MW:	33.9 kDa