

## Product datasheet for **RC219814L4V**

### **SIGLEC6 (NM\_198846) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

<b>Product Type:</b>	Lentiviral Particles
<b>Product Name:</b>	SIGLEC6 (NM_198846) Human Tagged ORF Clone Lentiviral Particle
<b>Symbol:</b>	SIGLEC6
<b>Synonyms:</b>	CD33L; CD33L1; CD33L2; CD327; CDW327; OBBP1
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Vector:</b>	pLenti-C-mGFP-P2A-Puro (PS100093)
<b>Tag:</b>	mGFP
<b>ACCN:</b>	NM_198846
<b>ORF Size:</b>	1059 bp
<b>ORF Nucleotide Sequence:</b>	The ORF insert of this clone is exactly the same as(RC219814).
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>RefSeq:</b>	<a href="#">NM_198846.1</a>
<b>RefSeq Size:</b>	2536 bp
<b>RefSeq ORF:</b>	1062 bp
<b>Locus ID:</b>	946
<b>UniProt ID:</b>	<a href="#">O43699</a>
<b>Cytogenetics:</b>	19q13.41
<b>Protein Families:</b>	Druggable Genome, Secreted Protein, Transmembrane
<b>MW:</b>	36.1 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the SIGLEC (sialic acid binding immunoglobulin-like lectin) family of proteins. The encoded transmembrane receptor binds sialyl-TN glycans and leptin. Placental expression of the encoded protein is upregulated in preeclampsia. [provided by RefSeq, Jul 2016]