

Product datasheet for RC219852L4

PAX3 (NM_181458) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: PAX3 (NM_181458) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: PAX3

Synonyms: CDHS; HUP2; WS1; WS3

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC219852).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





ACCN: NM_181458

ORF Size: 1449 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

PAX3 (NM_181458) Human Tagged Lenti ORF Clone - RC219852L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 181458.1</u>

RefSeq Size: 3359 bp
RefSeq ORF: 1455 bp
Locus ID: 5077
UniProt ID: P23760

 UniProt ID:
 P23760

 Cytogenetics:
 2q36.1

Protein Families: Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS,

Transcription Factors

MW: 53.3 kDa

Gene Summary: This gene is a member of the paired box (PAX) family of transcription factors. Members of the


PAX family typically contain a paired box domain and a paired-type homeodomain. These genes play critical roles during fetal development. Mutations in paired box gene 3 are associated with Waardenburg syndrome, craniofacial-deafness-hand syndrome, and alveolar rhabdomyosarcoma. The translocation t(2;13)(q35;q14), which represents a fusion between PAX3 and the forkhead gene, is a frequent finding in alveolar rhabdomyosarcoma. Alternative

splicing results in transcripts encoding isoforms with different C-termini. [provided by RefSeq,

Jul 2008]



Product images:



Circular map for RC219852L4