

Product datasheet for **RC219856L4V**

ZNF512 (NM_032434) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF512 (NM_032434) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF512
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032434
ORF Size:	1701 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219856).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032434.2
RefSeq Size:	3460 bp
RefSeq ORF:	1704 bp
Locus ID:	84450
UniProt ID:	Q96ME7
Cytogenetics:	2p23.3
Domains:	zf-C2H2
Protein Families:	Transcription Factors
MW:	64.5 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein containing four putative zinc finger motifs. Zinc finger motifs may bind to proteins or nucleic acids. Zinc finger-containing proteins are involved in a variety of processes, including regulation of transcription. Alternative splicing results in multiple transcript variants for this gene. [provided by RefSeq, Sep 2012]