

Product datasheet for **RC219858L4V**

NEU2 (NM_005383) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NEU2 (NM_005383) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NEU2
Synonyms:	SIAL2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005383
ORF Size:	1140 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219858).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005383.1 , NP_005374.1
RefSeq Size:	1143 bp
RefSeq ORF:	1143 bp
Locus ID:	4759
UniProt ID:	Q9Y3R4
Cytogenetics:	2q37.1
Protein Pathways:	Other glycan degradation, Sphingolipid metabolism
MW:	42.1 kDa



[View online »](#)

Gene Summary:

This gene belongs to a family of glycohydrolytic enzymes which remove sialic acid residues from glycoproteins and glycolipids. Expression studies in COS7 cells confirmed that this gene encodes a functional sialidase. Its cytosolic localization was demonstrated by cell fractionation experiments. [provided by RefSeq, Jul 2008]