

Product datasheet for **RC219876A1V**

Human KChIP2 (KCIP2) (NM_173191) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human KChIP2 (KCIP2) (NM_173191) AAV Particle
Tag:	Myc-DDK
Symbol:	KChIP2
Synonyms:	KCHIP2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC219876 representing NM_173191 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCGGGGCCAGGGCCGCAAGGAGAGTTTGTCCGATCCCCGAGACCTGGACGGCTCCTACGACCAGCTCA
CGGGCCACCCTCCAGGGCCCACTAAAAAGCGCTGAAGCAGCGATTCTCAAGCTGCTGCCGTGCTGCGG
GCCCCAAGCCCTGCCCTCAGTCAGTGAAACATTAGCCGCCCCAGCCTCCCTCCGCCCCACAGACCCCGC
CTGCTGGACCCAGACAGCGTGGACGATGAATTTGAATTGTCCACCGTGTGTACCGGCCCTGAGGGTCTGG
AGCAGCTGCAGGAGCAAACAAATTCACGCGCAAGGAGTTGCAGGTCTGTACCGGGGCTTCAAGAACGA
ATGTCCCAGCGGAATTGTCAATGAGGAGAACTTCAAGCAGATTTACTCCAGTTCTTTCTCAAGGAGAC
TCCAGCACCTATGCCACTTTTCTTTCAATGCCTTTGACACCAACCATGATGGCTCGGTCAAGTTTGAAG
ACTTTGTGGCTGGTTTGTCCGTGATTCTTCGGGGAAGTGTAGATGACAGGCTTAATTGGGCCTTCAACCT
GTATGACCTTAACAAGGACGGCTGCATCACCAGGAGGAAATGCTTGACATCATGAAGTCCATCTATGAC
ATGATGGGCAAGTACACGTACCCTGCACTCCGGGAGGAGGCCCAAGGGAACACGTGGAGAGCTTCTTCC
AGAAGATGGACAGAAACAAGGATGGTGTGGTGACCATTGAGGAATTCATTGAGTCTTGTCAAAAGGATGA
GAACATCATGAGGTCCATGCAGCTCTTTGACAATGTCATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC219876 representing NM_173191 Red=Cloning site Green=Tags(s) MRGQGRKESLSDSRDLGSDYDQLTGHPGPTKKALKQRFLKLLPCCGPQALPSVSETLAAPASLRPHRPR LLDPDSVDDEFELSTVCHRPEGLEQLQEQTKFTRKELQVLYRGFKNECPSGIVNEENFKQIYSQFFPQGD SSTYATFLFNAFDTNHDGSVSFEDFVAGLSVILRGTVDDRLNWFNLYDLNKDGCITKEEMLDIMKSIYD MMGKYTPALREEAPREHVESFFQKMDRNDGVTIEEFIESCQKDENIMRSMQLFDNVI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_173191
ORF Size:	810 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_173191.2</u>
RefSeq Size:	2563 bp
RefSeq ORF:	813 bp
Locus ID:	30819
UniProt ID:	<u>Q9NS61</u>
Cytogenetics:	10q24.32
MW:	30.7 kDa