

Product datasheet for **RC219880A1V**

Human beta 1 Sodium Potassium ATPase (ATP1B1) (NM_001001787) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human beta 1 Sodium Potassium ATPase (ATP1B1) (NM_001001787) AAV Particle
Tag:	Myc-DDK
Symbol:	beta 1 Sodium Potassium ATPase
Synonyms:	ATP1B; MGC1798
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC219880 representing NM_001001787 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCCGCGGAAAGCCAAGGAGGAGGGCAGCTGGAAGAAATTCATCTGGAACCTCAGAGAAGAAGGAGT
TTCTGGGCAGGACCGGTGGCAGTTGGTTTAAGATCCTTCTATTCTACGTAATATTTATGGCTGCCTGGC
TGGCATCTTCATCGGAACCATCCAAGTGATGCTGCTCACCATCAGTGAATTTAAGCCACATATCAGGAC
CGAGTGGCCCCGCCAGGATTAACACAGATTCTCAGATCCAGAAGACTGAAATTTCTTTCTGTCCTAATG
ATCCCAAGAGCTATGAGGCATATGTAAGTGAACATAGTTAGGTTCTGGAAAAGTACAAAGATTGAGCCCA
GAGGGATGACATGATTTTTGAAGATTGTGGCGATGTGCCAGTGAACCGAAAGAACGAGGAGACTTTAAT
CATGAACGAGGAGAGCGAAAGGTCTGCAGATTCAAGCTTGAATGGCTGGGAAATGCTCTGGATTAATG
ATGAACTTATGGCTACAAAGAGGGCAAACCGTGCATTATTATAAAGCTCAACCGAGTTCTAGGCTTCAA
ACCTAAGCCTCCCAAGAATGAGTCCTTGGAGACTTACCCAGTGATGAAGTATAACCCAAATGTCCTTCCC
GTTCACTGCACTGGCAAGCGAGATGAAGATAAGGATAAAGTTGGAAATGTGGAGTATTTGGACTGGGCA
ACTCCCTGGTTTTCTCTGCAGTATTATCCGTACTATGGCAACTCTGCAGCCCAATACCTGCAGCC
CCTGCTGGCCGTACAGTTCACCAATCTTACCATGGACACTGAAATTCGCATAGAGTGTAAAGCGTACGGT
GAGAACATTGGGTACAGTGAGAAAGACCGTTTTTCAGGGACGTTTTGATGTAAAAATTAATTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC219880 representing NM_001001787 Red=Cloning site Green=Tags(s) MARGKAKEEGSWKKFIWNSEKKEFLGRTGGSWFKILLFYVIFYGCLAGIFIGTIQVMLLTISEFKPTYQD RVAPPGLTQIPQIQKTEISFRNPKSYEAYVLNIVRFLEKYKDSAQRDDMIFEDCGDVPSEPKERGDFN HERGERKVCRFKLEWLGNCGLNDETYGYKEGKPCIIIKLNRVLGFKPKPPKNESLETYPVMKYNPVLP VQCTGKRDEDDKDKVGNVEYFGLGNSPGFPLQYYPPYGKLLQPKYLQPLLAVQFTNLTMDTEIRIECKAYG ENIGYSEKDRFQGRFDVKIKF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001001787
ORF Size:	903 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001001787.1, NP_001001787.1</u>
RefSeq Size:	1568 bp
RefSeq ORF:	905 bp
Locus ID:	481
Cytogenetics:	1q24.2
MW:	34.7 kDa