

Product datasheet for **RC219881L4V**

UQCRHL (NM_001089591) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	UQCRHL (NM_001089591) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UQCRHL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001089591
ORF Size:	273 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC219881).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001089591.1 , NP_001083060.1
RefSeq Size:	538 bp
RefSeq ORF:	276 bp
Locus ID:	440567
UniProt ID:	A0A096LP55
Cytogenetics:	1p36.21
Protein Pathways:	Alzheimer's disease, Cardiac muscle contraction, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease
MW:	10.6 kDa



[View online »](#)

Gene Summary:

This gene has characteristics of a pseudogene derived from the UQCRH gene. However, there is still an open reading frame that could produce a protein of the same or nearly the same size as that of the UQCRH gene, so this gene is being called protein-coding for now. [provided by RefSeq, Jul 2008]